Gravimetric Analysis:



Weighing the sample to be analyzed.	Dissolving this sample in water.	Adding a suitable chemical to form a precipitate.	Filtering to collect the precipitate	Repeated drying and weighing until a constant mass of precipitate is obtained.
---	--	---	--	--

Common Mistakes:

- Precipitate is not dry when you take the final mass.
 - Results in the appearance of more precipitate than was actually produced because some mass is water.
 - Percent yield would be higher than it should be.

Common Applications:

• Mixtures of solids—determining the amount of a particular ion in a solution

Important to Remember:

• All sodium, nitrate, ammonium, and potassium compounds are soluble. Net ionic equations would not include these ions.



Common Mistakes:

- Overfilling the volumetric flask
 - o Results in a dilute solution
- Not using distilled water.
 - o Other ions could affect the experiment for which the solution is used
- Not using a volumetric flask (beaker or Erlenmeyer instead)
 - o Loss of precision in concentration of prepared solution

Common Applications:

• Making solutions to dissolve substances for analysis, particularly in titrations.

Important to Remember:

• Molarity = moles solute/L of solution



Common Mistakes:

- Overshooting the titration (too dark of a color at the end)
 - Results in the concentration of the unknown solution in the flask appearing to be higher than it actually is, since too much titrant has been added.
- Not using indicator.
 - No perceivable endpoint.
- Using incorrect indicator.
 - pH at the equivalence point should be approximately equal to the pKa of the indicator.

Common Applications:

- Solving for the concentration of an unknown substance (analyte).
- Acid/Base, Redox

- Molarity = moles solute/L of solution
- Analyte: substance in flask
- Titrant: substance in buret
- Standard solution: solution of known concentration, usually goes into the buret.
- $M_1V_1 = M_2V_2$ is helpful for solving for the concentration of the analyte solution at the equivalence point.
- Endpoint: point in titration where flask solution changes color
- Equivalence point: point in the titration where the moles of acid are equal to the moles of base

Analyzing Concentration of Solutions Using Beer's Law):



Step 1: Pick the wavelength for the solution where absorbance is highest (for solute).



Step 2: Measure absorbance for different concentrations at that wavelength. Graph the results.



A = εbc

Absorbance = (molar absorptivity)(cuvette pathway length)(concentration)

Common Mistakes:

- Absorbance is lower than it should be (point falls below the line)
 - Cuvette was cleaned with distilled water and then immediately filled with solution, creating a more dilute solution
 - Too little solute in the prepared solution
- Absorbance is higher than it should be (point falls above the line)
 - \circ $\;$ Cuvette is dirty with fingerprints/dust, etc.
 - Too much solute in the prepared solution
 - o Contamination with a more concentrated solution
 - Used a cuvette with a longer path for one data point
- Did not use the correct wavelength of maximum absorbance for the solute.
 - Absorbances could be too low especially for dilute solutions
- Overfilled the cuvette
 - o Should not have an impact on data
- Picked a wavelength where it is high absorbance for the solvent
 - Won't be able to distinguish absorbance due to solvent vs. solute

Common Applications:

- Determining the concentration of a solution of unknown concentration using solutions of known concentration
- Kinetics reactions (like bleach + blue food dye)

- Before using, you need to calibrate the spectrophotometer with a blank of just solvent (in order to account for any absorbance due to solvent and cuvette itself)
- Molarity = moles solute/L of solution
- Absorbance is the amount of light the solution absorbs at a specific wavelength
- Molar absorptivity (1/M*cm) describes how intensely a sample absorbs light at a specific wavelength (constant unique to the substance at a specific wavelength)
- Path length of sample is the length of the cuvette where the light will travel (cm)
- Concentration is molarity

Chromatography





Common Mistakes:

- Solvent reaches the top of the paper strip.
 - Rf values cannot be calculated as we do not know how far the solvent would have traveled had their been more paper.
- No major difference in polarity between paper and solvent
 - Substances cannot be adequately separated
- No major differences in polarity of components of mixture
 - o Substances cannot be adequately separated

Common Applications:

• Determining the components of a mixture

- Paper is usually relatively nonpolar in comparison to the solvent.
- The substance that travels further up the paper is more attracted to the solvent.
- The substance that travels the least is most attracted to the paper.
- If multiple trials are run, compare Rf values, not relative heights.
- Polar substances tend to lack symmetry, have polar bonds, and have lone pairs on the central atom. They are most soluble in other polar substances.
- Nonpolar substances tend to be symmetrical, have identical bonds, and have no lone pairs on the central atom. They are most soluble in other nonpolar substances.

Fractional Distillation



Common Applications:

• Separating components in a solution/mixture based on differences in boiling point

- The substance with the lower boiling point has a greater vapor pressure and weaker intermolecular forces
- The substance with the higher boiling point has a lower vapor pressure and stronger intermolecular forces
- The temperature of the solution will remain constant while a component is boiling off

Coffee Cup Calorimetry



Common Mistakes:

• The final temperature is the highest (for exothermic) or lowest (for endothermic) temperature recorded during the reaction/process

Applications:

• Solving for the specific heat of a metal or the heat of reaction

- Endothermic processes have a drop in temperature.
- Exothermic processes have an increase in temperature.
- The water is not part of the system. It is part of the surroundings.
- $q = mC\Delta T$
 - q = heat in Joules or calories
 - o m = mass of entire solution OR object, grams or kilograms
 - \circ C = specific heat capacity, J/g°C (or a variation of the above)
 - $_{\circ}$ $\Delta T = T_{\text{final}} T_{\text{initial}}$