FLOWFLEX
HOME TEST

HOME TEST

FLOWFLEX
COVID-19 Antigen Home Test

HOME TEST

BEST PRACTICES

## **Tips for Flowflex:**

- This test takes 15 minutes to complete.
- Keep the box to use as a test stand. You will need to use it for the test.
- Use the easy to follow pictures and instructions inside the box on how to do the test.
- For those between the ages of 2 and 13, an adult may need to assist you in the test.
- Folow this QR Code to watch a video on how to do this test:



## **Report Your Test Results:**

 Report your results on Primary. Health! This information goes to your local Department of Public Health and School. It helps them keep track of cases in the community! Your information is safe and protected under this secure platform.

**How To Report:** 



**Already Registered:** 



Register & Report:

